定量NMR (qNMR)

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
図 1: General procedure of a qNMR experiment using an internal standard

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
図 2: Example of a ^1H qNMR spectrum with internal standard calibration.

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P_{Sample}	Purity of the sample as mass fraction
P_{CRM}	Purity of the CRM as mass fraction
I_{Analyte}	Integral of the analyte signal
I_{CRM}	Integral of the CRM signal
N_{Analyte}	Number of analyte protons (phosphorus nuclei, fluorine nuclei)
N_{CRM}	Number of CRM protons (phosphorus nuclei, fluorine nuclei)
M_{Analyte}	Molecular mass of the analyte
M_{CRM}	Molecular mass of the CRM
m_{Sample}	Mass of sample
m_{CRM}	Mass of CRM

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
図 3: Overview on typical contributions to the relative standard uncertainty (relative squared contributions are given) for the quantification of tris(2-chloroethyl) phosphate using phosphonoacetic acid as internal standard in a ^{31}P qNMR measurement.

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
Standards for ^1H quantitative NMR, TraceCERT®

Standards for ^{19}F quantitative NMR, TraceCERT®

Standards for ^{31}P quantitative NMR, TraceCERT®


Figure 4: Example spectra for all TraceCERT™ CRMs for qNMR in organic solution and / or D_2O respectively

TraceCERT™ Certified Reference Materials for qNMR

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TraceCERT®
Traceable Certified Reference Materials

SIGMA-ALDRICH®

Certificate

This certificate is designed in accordance with ISO Guide 31 [1].

Product name: **2,4-Dichlorobenzotrifluoride**

Product no.: **53396**

Lot no.: **BCBT0752**

Formula: **C₇H₃F₃Cl₂**

Molecular mass: **215.00 g/mol**

Traceability [2]: **NMIJ [3] CRM 4601-a**
(3,5-Bis(trifluoromethyl)benzoic acid)

Certificate issue date: **February 08, 2017**

Certificate version: **01**

Expiry date: **JAN 2019**

Traceability Statement

Expiry Date

Value (g/g)

Expanded Uncertainty (g/g)

2,4-Dichlorobenzotrifluoride

99.51 %

0.26 %

Expanded Uncertainty (g/g)

Minimum sample size: There is no recommended minimum sample weight. The substance is liquid at room-temperature and therefore homogeneous under these conditions. Potential inhomogeneity is covered by the expanded measurement uncertainty.

Intended use: Use this certified reference material (CRM) as internal standard for quantitative ^{19}F -NMR measurements only. The content of this CRM was determined by ^{19}F -qNMR spectroscopy.

Storage and handling: The CRM should be stored in the original bottle at room-temperature (20-25°C). After use the bottle should be tightly closed and protected from excessive moisture and light. Storage under Argon is recommended.

CRM operations: *A. Rück*

Certification body: *Dr. K.D. Schmidt*

Certification page 1 of 4

Sigma-Aldrich Production GmbH, Industriestrasse 25, 9471 Buchs/Schweizland, Tel: +41 81-756-2011, Fax: +41 81-756-5449

図 5: Example of a Certificate 94681 Methyl 3,5-dinitrobenzoate (first page)

References

1. F. Malz, H. Jancke, J. Pharm. Biomed. Anal. 38 (2005), 813–23.
2. S.K. Bharti, R. Roy, Quantitative ^1H NMR spectroscopy, Trends Anal. Chem. 35(2012) 5–26.
3. M. Weber, C. Hellriegel, A. Rueck, R. Sauermoser, J. Wuethrich, Using high-performance quantitative NMR (HP-qNMR®) for certifying traceable and highly accurate purity values of organic reference materials with uncertainties <0.1 %, Accred. Qual. Assur. 18 (2013), 91–98.
4. T. Saito, T. Ihara, M. Koike, S. Kinugasa, Y. Fujimine, K. Nose, T. Hirai, A new traceability scheme for the development of international system-traceable persistent organic pollutant reference materials by quantitative nuclear magnetic resonance, Accred. Qual. Assur. 14 (2009) 79–86.
5. T. Ihara, T. Saito, N. Sugimoto, Expansion of organic reference materials for the analysis of hazardous substances in food and the environment, Synthesiol. Engl. Ed. 2 (2009) 12–22.
6. Eurachem/CITAC Guide, 2nd Ed. (2000), "Quantifying uncertainty in analytical measurement."
7. M. Weber, C. Hellriegel, A. Rueck, J. Wuethrich, P. Jenks, Using high-performance ^1H NMR (HP-qNMR®) for the certification of organic reference materials under accreditation guidelines - Describing the overall process with focus on homogeneity and stability assessment, JPBA 93 (2014), 102–110.
8. M. Weber, C. Hellriegel, A. Rueck, J. Wuethrich, P. Jenks, M. Obkircher, Analytical and Bioanalytical Chemistry (2015), 407(11), 3115–3123.

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本紙記載の製品は試験・研究用です。ヒト、動物への治療、もしくは診断目的として使用しないようご注意ください。掲載価格は希望販売価格(税別)です。実際の価格は弊社製品取扱販売店へご確認ください。なお、品目、製品情報、価格等は予告なく変更される場合がございます。予めご了承ください。記載内容は2021年5月時点の情報です。Merck, the vibrant M, and Supelco are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. ©2021 Merck KGaA, Darmstadt, Germany. All rights reserved.