

機能性高分子材料の 精密重合を可能にする

RAFT剤

トリチオカルボナート
ジチオカルバマート
ジチオベンゾアート
切り替え可能なRAFT剤
マクロRAFT剤

トリチオカルボナート

構造式	化合物名	適用可能なモノマー	カタログ番号	Ref.
	Cyanomethyl dodecyl trithiocarbonate, 98%	styrenes, acrylates, acrylamides	723029-1G 723029-5G	1, 35-37 34 (化合物171)
	4-Cyano-4-[(dodecylsulfanylthiocarbonyl)sulfanyl]pentanoic acid, 97%	methacrylates, methacrylamides, styrenes	723274-1G 723274-5G	2, 14, 35-37, 41 34 (化合物98)
	2-(Dodecylthiocarbonothioylthio)propionic acid, 97%	methacrylates, methacrylamides	749133-1G 749133-5G	3, 4, 14, 37 34 (化合物155)
	2-Cyano-2-propyl dodecyl trithiocarbonate, 97%	methacrylates, methacrylamides, styrenes	723037-1G 723037-5G	5, 6, 35-37 34 (化合物95)
	2-(Dodecylthiocarbonothioylthio)-2-methylpropionic acid 3-azido-1-propanol ester, 98% アジド基は、様々なアルキン官能基化生体分子とのコンジュゲートに使用可能です。	styrenes, acrylates, acrylamides	741698-1G 741698-5G	7, 37 34 (化合物137)
	4-Cyano-4-[(dodecylsulfanylthiocarbonyl)sulfanyl]pentanol このアルコール官能基化RAFT剤は、標準的なカップリング反応や開環重合など、重合前および重合後に様々な官能基化を行うことができ、十分に制御されたブロック共重合体を得ることができます。	methacrylates, methacrylamides, styrenes	760110-1G 760110-5G	8, 37
	2-(Dodecylthiocarbonothioylthio)-2-methylpropionic acid, 98%	styrenes, acrylates, acrylamides	723010-1G 723010-5G	9, 36, 37, 41 34 (化合物123)
	Methyl 2-(dodecylthiocarbonothioylthio)-2-methylpropionate, 97%	styrenes, acrylates, acrylamides	740497-5G	10, 37

シグマ アルドリッチのバルクサービス

カタログ掲載最大容量の
20倍よりバルクのご相談が可能です。

バルク対応お問い合わせフォーム www.sigmaaldrich.com/chem-bulk-jp
からお気軽にお問い合わせください。



構造式	化合物名	適用可能なモノマー	カタログ番号	Ref.
	2-(Dodecylthiocarbonothioylthio)-2-methylpropionic acid <i>N</i> -hydroxysuccinimide ester, 98% NHSエステル末端基は、様々な生体分子との結合に使用されます。	styrenes, acrylates, acrylamides	741035-1G	11, 37
	Pentaerythritol tetrakis[2-(dodecylthiocarbonothioylthio)-2-methylpropionate], 97%	styrenes, acrylates, acrylamides	763551-1G	47, 48
	<i>S,S</i> -Dibenzyl trithiocarbonate, 97%	styrenes, methacrylates, acrylamides	746304-1G 746304-5G	12, 13, 49 34(化合物88)
	2-(2-Carboxyethylsulfanylthiocarbonylsulfanyl) propionic acid 水溶性かつ有機溶媒可溶性の両親媒性RAFT剤	acrylates, acrylamides, styrenes	900152-1G	15 34(化合物157)
	4-((((2-Carboxyethyl)thio)carbonothioyl)thio)-4-cyanopentanoic acid, 95% 水溶性かつ有機溶媒可溶性の両親媒性RAFT剤	acrylates, acrylamides, styrenes, methacrylates, methacrylamides	900161-1G	15
	Bis(dodecylsulfanylthiocarbonyl) disulfide, 95+% RAFT剤合成のための前駆体化合物		723126-5G	16, 35



ジチオカルバマート

構造式	化合物名	適用可能なモノマー	カタログ番号	Ref.
	Cyanomethyl (3,5-Dimethyl-1 <i>H</i> -pyrazole)-carbodithioate, 95% ピラゾール系RAFT剤は、有機溶媒に可溶であり、無臭です。電子的に全く異なる多様なモノマーのグループを使用した多機能性ポリマーの合成が可能です。	acrylates, acrylamides, styrenes	900150-1G	17, 18, 50
	2-Cyanobutan-2-yl 4-chloro-3,5-dimethyl-1 <i>H</i> -pyrazole-1-carbodithioate, 95% 「多用途RAFT剤」と称され、全てのモノマーの重合を行うことができます。多分散度が重要な用途での使用が推奨されます。	acrylates, acrylamides, styrenes, methacrylates, methacrylamides, vinyl esters, vinyl amides	900158-1G	17, 18, 51
	Cyanomethyl methyl(phenyl) carbamodithioate, 98%	vinyl esters, vinyl amides	723002-1G	20, 21, 36-38 34(化合物258)











ジチオベンゾアート

構造式	化合物名	適用可能なモノマー	カタログ番号	Ref.
	1-(Methoxycarbonyl)ethyl benzodithioate, 97+%	methacrylates, methacrylamides	751138-1G	37, 52
	2-Cyano-2-propyl benzodithioate, >97%	methacrylates, methacrylamides	722987-1G 722987-5G	19, 22, 23, 35-37 34(化合物12)
	4-Cyano-4-(phenylcarbonothioylthio) pentanoic acid	methacrylates, methacrylamides	722995-1G 722995-5G	24, 25, 36, 37 34(化合物18)
	4-Cyano-4-(phenylcarbonothioylthio) pentanoic acid <i>N</i> -succinimidyl ester	methacrylates, methacrylamides	758353-1G	26-28, 37
	Ethyl 2-(phenylcarbonothioylthio)-2-phenylacetate, 98% メチルエステル体よりも高い温度安定性を示します。	methacrylates, methacrylamides	773506-1G	29, 30
	Bis(thiobenzoyl) disulfide, >95% RAFT剤合成のための前駆体化合物		723118-5G	16, 35

切り替え可能な/汎用 (switchable/universal) RAFT 剤

構造式	化合物名	適用可能なモノマー	カタログ番号	Ref.
	2-Cyanopropan-2-yl <i>N</i> -methyl- <i>N</i> -(pyridin-4-yl) carbamodithioate, 97%	脱プロトン化: vinyl esters, vinyl amides プロトン化: styrenes, acrylates, methacrylate	736236-1G	31, 39, 53, 54
	<i>N,N'</i> -Dimethyl <i>N,N'</i> -di(4-pyridinyl)thiuram disulfide 切り替え可能なRAFT剤合成用前駆体化合物		735973-5G	39, 54

マクロRAFT剤

構造式	化合物名	適用可能なモノマー	カタログ番号	Ref.
	Poly(ethylene glycol) methyl ether (4-cyano-4-pentanoate dodecyl trithiocarbonate), average Mn 5,400	methacrylates, methacrylamides, styrenes	751626-1G	37
	Poly(ethylene glycol) methyl ether 4-cyano-4-[(dodecylsulfanylthiocarbonyl)sulfanyl] pentanoate, average Mn 10,000	methacrylates, methacrylamides, styrenes	753033-1G	34 (化合物132)
	Poly(ethylene glycol) bis[2-(dodecylthiocarbonothioylthio)-2-methylpropionate], average Mn 10,800	styrene, acrylate, acrylamide	753025-1G	37, 55
	Poly(ethylene glycol) methyl ether 2-(dodecylthiocarbonothioylthio)-2-methylpropionate, average Mn 6,000	styrene, acrylate, acrylamide	736325-1G	37, 44, 46, 55
	Poly(ethylene glycol) 4-cyano-4-(phenylcarbonothioylthio)pentanoate, average Mn 2,000	methacrylates, methacrylamides	764914-1G	45 34 (化合物19)
	Poly(ethylene glycol) 4-cyano-4-(phenylcarbonothioylthio)pentanoate, average Mn 10,000		764930-1G	
	Poly(<i>N,N</i> -dimethylacrylamide), DDMAT terminated, average Mn 10,000, PDI ≤ 1.1 末端のDDMAT構造を用いて、バイオマテリアルやリソグラフィー用途で重要なブロック共重合体を作ることができます。	styrene, acrylate, acrylamide	773638-1G	42
	Polystyrene, DDMAT terminated, average Mn 5,000, PDI < 1.1	styrene, acrylate, acrylamide	772577-1G	43
	Polystyrene- <i>block</i> -poly(acrylic acid), DDMAT terminated, PS:PAA 3,000:5,000, PDI ≤ 1.1 トリブロック共重合体の生成に使用可能なRAFT構造末端基を有しており、ライゲーシンのための修飾も可能です。	styrene, acrylate, acrylamide	776351-500MG	40
	Bis-MPA-RAFT dendrimer, trimethylol propane core, generation 1 このデンドリマーには6つのRAFT構造末端基が含まれており、メタクリル系モノマーの制御されたラジカル重合により6本鎖のスターポリマーを合成することができます。	acrylate	911550-250MG	32, 33

精密ラジカル重合ガイドブック

カタログPDFのダウンロードはこちら
<https://bit.ly/crp-guide>



- 原子移動ラジカル重合 (ATRP)
- 可逆的付加-開裂連鎖移動重合 (RAFT)
- ニトロキンドを介した重合 (NMP)
- モノマー製品一覧



References

1. Dual Stimuli-Responsive Self-Assembly Behavior of a Tailor-Made ABC-Type Amphiphilic Tri-Block Copolymer. Pal S., et al., *J. Polym. Sci.* **2020**, *58*, 843–851.
2. Advances in RAFT polymerization: the synthesis of polymers with defined end-groups. G. Moad, et al., *Polymer* **2005**, *46*, 8458–8468.
3. Block copolymers of poly(2-oxazoline)s and poly(meth) acrylates: a crossover between cationic ring-opening polymerization (CROP) and reversible addition-fragmentation chain transfer (RAFT). Krieg, A., et al. *ACS Macro Letters* **2012**, *1*, 776–779.
4. Carbon fibers prepared from tailored reversible-addition-fragmentation transfer copolymerization-derived poly(acrylonitrile)-co-poly(methylmethacrylate). Spori J.M., et al., *J. Polym. Sci., Part A: Polym. Chem.* **2014**, *52*, 1322–1333.
5. A novel solid state photocatalyst for living radical polymerization under UV irradiation. Fu, Q., McKenzie, T., Ren, J. et al., *Sci Rep.* **2016**, *6*, 20779.
6. Photo-induced reversible addition-fragmentation chain transfer (RAFT) polymerization of acrylonitrile at ambient temperature: A simple system to obtain high-molecular-weight polyacrylonitrile. Li J., et al., *React. Funct. Polym.* **2017**, *113*, 1–5.
7. Versatile Pathway to Functional Telechelics via RAFT Polymerization and Click Chemistry. Gondi S.R., et al., *Macromolecules* **2007**, *40*, 474–481.
8. One-step synthesis of block copolymers using a hydroxyl-functionalized trithiocarbonate RAFT agent as a dual initiator for RAFT polymerization and ROP. Kang, H. U., et al., *J. Polym. Sci. Part A: Polym. Chem.* **2013**, *51*, 774–779.
9. Functional Polymers from Novel Carboxyl-Terminated Trithiocarbonates as Highly Efficient RAFT Agents. Lai, J.T., et al., *Macromolecules* **2002**, *35*, 6754–6756.
10. Synthesis of PMMA Microparticles with a Narrow Size Distribution by Photoinitiated RAFT Dispersion Polymerization with a Macromonomer as the Stabilizer. Tan J., et al., *Macromolecules* **2014**, *47*, 6856–6866.
11. Construction of polymer coated core-shell magnetic mesoporous silica nanoparticles with triple responsive drug delivery. Hegazy M., et al., *Polym. Chem.* **2017**, *8*, 5852–5864.
12. RAFT Agent Design and Synthesis. Keddie, D. J., et al., *Macromolecules* **2012**, *45*, 5321–5342.
13. Dibenzy trithiocarbonate mediated reversible addition-fragmentation chain transfer polymerization of acrylonitrile. Liu X.-H., et al., *Inc. J. Polym. Sci. Part A: Polym. Chem.* **2006**, *44*, 490–498.
14. Selective Photoactivation of Trithiocarbonates Mediated by Metal Naphthalocyanines and Overcoming Activation Barriers Using Thermal Energy. Wu, Z., et al., *J. Am. Chem. Soc.* **2022**, *144*, 2, 995–1005.
15. Synthesis of Discrete Oligomers by Sequential PET-RAFT Single-Unit Monomer Insertion. Xu J. et al., *Angew. Chem. Int. Ed.* **2017**, *56*, 8376–8383.
16. A novel synthesis of functional dithioesters, dithiocarbamates, xanthates and trithiocarbonates. Tang S. H., et al., *Tetrahedron Lett.* **1999**, *40*, 2435–2438.
17. Dithiocarbamate RAFT agents with broad applicability - the 3,5-dimethyl-1H-pyrazole-1-carbodithioates. Gardiner J, et al., *Polym. Chem.* **2016**, *7*, 481–492.
18. 4-Halogeno-3,5-dimethyl-1H-pyrazole-1-carbodithioates: versatile reversible addition fragmentation chain transfer agents with broad applicability. Gardiner J. et al., *Polym. Int.* **2017**, *66*, 1438–1447.
19. The pH-responsive behaviour of poly(acrylic acid) in aqueous solution is dependent on molar mass. Swift T., et al., *Soft Matter* **2016**, *12*, 2542–2549.
20. Reversible Addition-Fragmentation Chain Transfer Polymerization of Vinyl Chloride. Abreu C.M., et al., *Macromolecules* **2012**, *45*, 2200–2208.
21. Synthesis of Defined Polymers by Reversible Addition—Fragmentation Chain Transfer: The RAFT Process. Rizzardo, E., et al., *ACS Symp. Ser.* **2000**, *768*, 278–96.
22. Reversible Addition Fragmentation Chain Transfer Polymerization of Methyl Methacrylate in the Presence of Lewis Acids: An Approach to Stereocontrolled Living Radical Polymerization. Chong Y.K., et al., *Macromolecules* **2007**, *40*, 9262–9271.
23. Fluorescent and chemico-fluorescent responsive polymers from dithiomaleimide and dibromomaleimide functional monomers. Robin M.P. and Reilly R.K. *Chem. Sci.* **2014**, *5*, 2717–2723.
24. Aqueous RAFT Polymerization of Acrylonitrile. Kopec M., et al., *Macromolecules* **2016**, *49*, 5877–5883.
25. One-pot synthesis of hyperbranched glycopolymers by RAFT polymerization. Wei Z., et al., *J. Polym. Sci. Part A: Polym. Chem.* **2012**, *50*, 2378–2388.
26. Polyethylene Building Blocks by Catalyzed Chain Growth and Efficient End Functionalization Strategies, Including Click Chemistry. Briquel, R., et al., *Angew. Chem. Int. Ed.*, **2008**, *47*, 9311–9313.
27. A switchable polymer brush system for antifouling and controlled detection. Demirci S., et al., *Chem. Commun.* **2017**, *53*, 3713–3716.
28. Dissipative Self-Assembly of Dynamic Multicompartmentalized Microsystems with Light-Responsive Behaviors. Cheng G., et al., *Chem* **2020**, *6*, 1160–1171.
29. Effect of the Backbone Tether on the Electrochemical Properties of Soluble Cyclopropenium Redox-Active Polymers. Montoto E.C., et al., *Macromolecules* **2018**, *51*, 3539–3546.
30. A facile random copolymer strategy to achieve highly conductive polymer gel electrolytes for electrochemical applications. Kim Y.M., et al., *J. Mater. Chem. C* **2019**, *7*, 161–169.
31. Universal (Switchable) RAFT Agents. Benaglia, M., et al., *J. Am. Chem. Soc.* **2009**, *131*, 6914–6915.
32. Dendritic architectures based on bis-MPA: functional polymeric scaffolds for application-driven research. Carlmark A., et al., *Chem. Soc. Rev.*, **2013**, *42*, 5858–5879.
33. Chemistry of multifunctional polymers based on bis-MPA and their cutting-edge applications. Gallego S.G., et al., *Prog. Polym. Sci.* **2015**, *48*, 85–110.
34. Living Radical Polymerization by the RAFT Process – A Third Update. Mode G., et al., *Aust. J. Chem.* **2012**, *65*, 985.
35. シグマ アルドリッチ 精密ラジカル重合ガイドブック;「高活性なモノマー」(MAM)のRAFT重合, p22
36. シグマ アルドリッチ 精密ラジカル重合ガイドブック; RAFTにおける典型的な重合方法, p25
37. シグマ アルドリッチ 精密ラジカル重合ガイドブック; RAFT剤とモノマーの適合表, p50–51
38. シグマ アルドリッチ 精密ラジカル重合ガイドブック;「低反応性モノマー」(LAM)のRAFT重合, p23
39. シグマ アルドリッチ 精密ラジカル重合ガイドブック; 汎用/切り替え可能なRAFT剤を用いる重合方法, p27
40. Multiple Morphologies of PAA-*b*-PSt Assemblies throughout RAFT Dispersion Polymerization of Styrene with PAA Macro-CTA. He W.-D., et al., *Macromolecules* **2011**, *44*, 9, 3358–3365.
41. Bilayer-folded lamellar mesophase induced by random polymer sequence. Shin, M., et al., *Nat. Commun.* **2022**, *13*, 2433.
42. Surfactant-Free RAFT Emulsion Polymerization Using Poly(*N,N*-dimethylacrylamide) Trithiocarbonate Macromolecular Chain Transfer Agents. Rieger J., et al., *Macromolecules* **2010**, *43*, 6302–6310.
43. A well-defined polyelectrolyte and its copolymers by reversible addition fragmentation chain transfer (RAFT) polymerization: synthesis and applications. Mumtaz M., et al., *RSC Adv.* **2015**, *5*, 98559–98565.
44. Adding a solvophilic comonomer to the polymerization-induced self-assembly of block copolymer and homopolymer: a cooperative strategy for preparing large compound vesicles. Huang C., et al., *RSC Adv.* **2017**, *7*, 46069–46081.
45. RAFT Aqueous Dispersion Polymerization Yields Poly(ethylene glycol)-Based Diblock Copolymer Nano-Objects with Predictable Single Phase Morphologies. Warren N. J., et al., *J. Am. Chem. Soc.* **2014**, *136*, 1023–1033.
46. In Situ Synthesis of Block Copolymer Nanoassemblies via Polymerization-Induced Self-Assembly in Poly(ethylene glycol). Gao C., et al., *Macromolecules* **2016**, *49*, 3789–3798.
47. PET-RAFT facilitated 3D printable resins with multifunctional RAFT agents. Shi, X., et al., *Mater. Chem. Front.* **2021**, *5*, 2271–2282.
48. Pressure-assisted strategy for the synthesis of vinyl pyrrolidone-based macro-star photoinitiators. A route to star block copolymers. Maksym, P., et al., *J. Polym. Sci.* **2020**, *58*, 1393–1399.
49. A Photoinduced Dual-Wavelength Approach for 3D Printing and Self-Healing of Thermosetting Materials. Zhang, Z., et al., *Angew. Chem. Int. Ed.* **2022**, *61*, e202114111.
50. Effect of Thiocarbonylthio Compounds on Visible-Light-Mediated 3D Printing. Zhang, Z., et al., *Macromolecules* **2021**, *54*, 3, 1170–1182.
51. Precise Control of Both Dispersity and Molecular Weight Distribution Shape by Polymer Blending. Whitfield, W., et al., *Angew. Chem. Int. Ed.* **2021**, *60*, 19383–19388.
52. Reversible addition-fragmentation chain transfer radical copolymerization of β -pinene and methyl acrylate. Wang, Y., et al., *Eur. Polym. J.* **2006**, *42*, 2695–2702.
53. Synthesis of CO₂-responsive gradient copolymers by switchable RAFT polymerization and their controlled self-assembly. Guo, X., et al., *Polym. Chem.* **2020**, *11*, 6794–6802.
54. High Rates and Substrate Selectivities in Water by Polyvinylimidazoles as Transaminase Enzyme Mimics with Hydrophobically Bound Pyridoxamine Derivatives as Coenzyme Mimics. Skouta, R., et al., *J. Am. Chem. Soc.* **2009**, *131*, 43, 15604–15605.
55. Molecular Design of Chemically Fueled Peptide–Polyelectrolyte Coacervate-Based Assemblies. Späth, F., et al., *J. Am. Chem. Soc.* **2021**, *143*, 12, 4782–4789.



本紙記載の製品は試験・研究用です。ヒト、動物への治療、もしくは診断目的として使用しないようご注意ください。掲載価格は希望販売価格(税別)です。実際の価格は弊社製品取扱販売店へご確認ください。なお、品目、製品情報、価格等は予告なく変更される場合がございます。予めご了承ください。記載内容は2025年3月時点の情報です。Merck, the vibrant M, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. ©2025 Merck KGaA, Darmstadt, Germany. All rights reserved.

シグマ アルドリッチ ジャパン

ライフサイエンス リサーチソリューションズ事業部

〒106-0041 東京都港区麻布台1-3-1 麻布台ヒルズ 森JPタワー 26階

製品の最新情報はこちら www.sigmaaldrich.com/JP/ja

製品に関するお問い合わせは、テクニカルサービスへ

E-mail: jpts@merckgroup.com Tel: 03-6756-8245

在庫照会・ご注文に関するお問い合わせは、カスタマーサービスへ

E-mail: sialjpcs@merckgroup.com Tel: 03-6756-8275

シグマ アルドリッチ ジャパン合同会社はメルクのグループ会社です。