

ケムビーズ (ChemBeads) 触媒コーティングしたガラスビーズ



スクリーニングを より早く・最適に

ケムビーズは、ミリグラム未満の固体を効率的かつ正確に供給する必要がある作業、特にハイスループット実験や並列合成用に作られた、触媒をコーティングした小さなガラスビーズです。ケムビーズ製品は通常のガラスビーズと同様の挙動をするため高い均一性や流動性を有し、凝集を防止することができます。

● 計量の簡素化

較正済みの専用スパチュラを使用して使用量を測定

● 時間の節約

ケムビーズはロボットプラットフォームで使用でき、1回のプログラミングで製品ライン全体を自動化

● 触媒の無駄を削減

ナノモルレベルの量を正確に添加可能



ケムビーズの詳細はこちら

<http://SigmaAldrich.com/ChemBeads>



ご注文情報

製品名	製品番号
Grubbs Catalyst® M204 ChemBeads	919764
Tris(dibenzylideneacetone)dipalladium(0) ChemBeads	919772
[1,1'-Bis(diphenylphosphino)ferrocene] dichloropalladium(II) ChemBeads	919780
Palladium(II) chloride ChemBeads	919853
Palladium(II) Acetate ChemBeads	924377
Bis(triphenylphosphine)palladium(II) dichloride ChemBeads	927759
XPhos Pd G2 ChemBeads	927775
XPhos Pd G3 ChemBeads	927783
Bis(di-tert-butyl(4-dimethylaminophenyl)phosphine) dichloropalladium(II) ChemBeads	927791
Zinc Acetate ChemBeads	927805
RuPhos Pd G3 ChemBeads	928283
Platinum(IV) oxide ChemBeads	928291
cataCXium® A Pd G3 ChemBeads	928305
tBuBrettPhos Pd G3 ChemBeads	928313
BrettPhos Pd G3 ChemBeads	928321
BrettPhos ChemBeads	928348
Xantphos ChemBeads	928356
XPhos ChemBeads	928364
Pd-PEPPSI™-IPent Catalyst ChemBeads	928399
tBuXPhos Pd G3 ChemBeads	931055
PEPPSI™-IPr catalyst ChemBeads	931063
Potassium hexachloroplatinate(IV) CHEMBEADS	932175
Hexaammineruthenium(III) chloride ChemBeads	932183
SPhos ChemBeads	932191
RuPhos ChemBeads	932205
XantPhos Pd G3 ChemBeads	932213
Grubbs Catalyst® M102 Chembeads	936685
Hoveyda-Grubbs Catalyst® M721 Chembeads	936723
tBuBrettPhos ChemBeads	936758
JosiPhos SL-J009-1 Pd G3 ChemBeads	936766
PPh3 Pd G2 ChemBeads	936774
P(oTol)3 Pd G2 ChemBeads	936782
meCgPPh Pd G3 ChemBeads	936790
APhos Pd G3 ChemBeads	936820
RockPhos Pd G3 Chembeads	936863
AdBrettPhos Pd G3 ChemBeads	936871
SPhos Pd G3 ChemBeads	936898

製品名	製品番号
XPhos Pd G4 ChemBeads	936928
SPhos Pd G4 ChemBeads	936936
RuPhos Pd G4 ChemBeads	936944
BrettPhos Pd G4 ChemBeads	936952
DavePhos Pd G3 Chembeads	936960
rac-BINAP Pd G3 ChemBeads	936979
DTBPF PD G3 CHEMBEADS	936987
PCy3 Pd G4 Chembeads	937002
XantPhos Pd G4 ChemBeads	937010
rac-BINAP Pd G4 Chembeads	937029
cataCXium® A Pd G3	937037
APhos Pd G4 ChemBeads	937088
P(t-Bu)3 Pd G4 ChemBeads	937096
GPhos ChemBeads	938750
Cesium carbonate Chembeads	938777
Sodium tert-butoxide Chembeads	938785
Potassium Phosphate Chembeads	938793
Potassium Fluoride Chembeads	938807
Potassium tert-butoxide Chembeads	938963
Sodium hydroxide Chembeads	938998
Cesium Fluoride Chembeads	939005
Copper(I) iodide Chembeads	939013
Sodium carbonate Chembeads	939021
Sodium phosphate Chembeads	939056
Lithium hydroxide Chembeads	939064

較正済みの専用スパチュラで計量も簡単に

製品名	製品番号
Calibrated ChemBeads Scoop	925470



メルクライフサイエンス公式 SNS、動画コンテンツをご覧ください。

本紙記載の製品は試験・研究用です。ヒト、動物への治療、もしくは診断目的として使用しないようご注意ください。掲載価格は希望販売価格(税別)です。実際の価格は弊社製品取扱販売店へご確認ください。なお、品目、製品情報、価格等は予告なく変更される場合がございます。予めご了承ください。記載内容は2025年3月時点の情報です。Merck, the vibrant M, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. ©2025 Merck KGaA, Darmstadt, Germany. All rights reserved.

シグマ アルドリッチ ジャパン

ライフサイエンス リサーチソリューションズ事業部

〒106-0041 東京都港区麻布台1-3-1 麻布台ヒルズ 森JPタワー 26階

製品の最新情報はこちら www.sigmaaldrich.com/JP/ja

製品に関するお問い合わせは、テクニカルサービスへ

E-mail: jpts@merckgroup.com Tel: 03-6756-8245

在庫照会・ご注文に関するお問い合わせは、カスタマーサービスへ

E-mail: sialjpcs@merckgroup.com Tel: 03-6756-8275

シグマ アルドリッチ ジャパン合同会社はメルクのグループ会社です。