



## Explore further with New Chemical Synthesis Products

We remain dedicated to supporting all your explorations by making the latest innovations in chemistry accessible in a bottle. We partner with the brightest minds in chemistry and have highlighted a few of our newest products below. Visit our website to view all of the products that will help you explore more and suggest any addition to our catalog that will help you even more.




### Micellar Catalysis

EvanPhos is a new biaryl phosphine ligand designed in the Lipshutz Lab for use in Suzuki-Miyaura cross couplings in conjunction with a Pd(0) source. Catalyst loadings as low as 0.1 mol% have been proven effective even when coupling densely functionalized partners. EvanPhos can be used in either organic solvent or under aqueous micellar conditions with the surfactant TPGS-750-M, making this an attractive technology for greener chemistry.



### Electrophilic Fluorination

Fluorination is an invaluable tool for medicinal chemistry as aryl fluorides are widely used in the pharmaceutical industry. This new Pd catalyst from the Ritter group enables direct C-H fluorination of unactivated arenes using either NFSI or Selectfluor as the fluoronium source.




### C-H Activation


This 2-pyridone ligand, which was developed by the laboratory of Jin-Quan Yu, binds to palladium and accelerates non-directed C-H functionalization. Developed for C-H olefination and carboxylation with arene as the limiting reagent, the Wang-Yu ligand enables diversification of drugs, synthetic intermediates, and other bioactive small molecules.

To view these and other new products, visit [SigmaAldrich.com/newchemistry](https://www.sigmaaldrich.com/newchemistry).


## Catalysis




905240



902284




904937




905534


## Chemical Biology




901222



903825




901183




901848


## Reagents




903132



900565




902918




905860


## Building Blocks




902616



902713



901783



901749

To view these and other new products, visit [SigmaAldrich.com/newchemistry](https://www.sigmaaldrich.com/newchemistry).

To place an order or receive technical assistance in Europe, please call Customer Service:  
 France: 0825 045 645 Spain: 901 516 645 Option 1  
 Germany: 069 86798021 Switzerland: 0848 645 645  
 Italy: 848 845 645 United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840  
 Or visit: [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices)  
 For Technical Service visit: [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice)

Merck KGaA  
 Frankfurter Strasse 250  
 64293 Darmstadt, Germany

[SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2019 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, Sigma-Aldrich, and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

